

Start M2 Effect

The all-rounder in the Start family




The Start M2 Effect is the all-rounder of the Start product family. This basic model has the same frame module as the more complex versions.

With its aluminium frame and a uniform caster fork for all wheel sizes, the Start M2 Effect makes a straightforward impression. Its outstanding features include a narrow folding size and durable upholstery.

Thanks to its modular design, the Start M2 Effect can be quickly and easily adapted to your individual requirements at any time.



- The side panel can be swung completely out of the way for easy transfers.



- The rear wheel adapter can be adjusted to a short/active or long/passive wheelbase.



- The aluminium crossbrace provides stability and good folding ability.

Options



• Start M2 Effect with z10-ce add-on drive



• Folding back tubes for easy and space-saving transportation



• Side panel with height-adjustable arm rest



• Angle-adjustable back can be easily adjusted without tools

Product Features

- Aluminium tube frame
- All main component groups are attached with screws
- Standardised caster forks for all wheel sizes
- Durable upholstery
- Narrow folding size with folding fixture
- Sturdy one-and-half fold crossbrace
- Knee lever wheel lock – standard



Basic Technical Data

Load capacity	125 kg
Weight	From 15 kg
Seat width	38 cm - 50.5 cm
Seat depth	40 cm, 43 cm, 46 cm
Seat height	42.5 cm, 47.5 cm
	Front 42 cm - 51 cm Rear 38 cm - 49 cm
Lower leg length	27 cm - 47 cm