

Start M4 XXL


The robust lightweight wheelchair



The Start M4 XXL is the robust model among the lightweight wheelchairs. With its sturdy double crossbrace and the reinforced rear wheels, it can carry up to 160 kilograms and withstand extreme stress and strain.

The optional back stabiliser rounds out the wheelchair's features without compromising its foldability.

The rear wheel adapter of the Start M4 XXL allows the centre of gravity to be adjusted to meet different user needs.



- The continuous aluminium foot plate ensures greater stability and meets maximum space requirements.



- Enhanced safety when driving downhill - and the optional drum brake system for the attendant.



- The rear wheel adapter permits various wheelbase settings.

Options



• Start M4 with z10 push assist



• Optional swinging anti-tipper allows to overcome low obstacles without deactivation



• Back stabiliser improves stability and can be swung out of the way for folding



• Double crossbrace permits a maximum load capacity of 160 kg

Product Features

- Aluminium tube frame
- All main component groups are attached with screws
- Reinforced rear wheels and back tubes
- Optional stabiliser bar
- Standardised caster forks for all wheel sizes
- Durable upholstery
- Narrow folding size with folding fixture
- Sturdy double crossbrace
- Optimised knee lever wheel lock



Basic Technical Data

Load capacity	160 kg
Weight	From 17.5 kg
Seat width	43 cm - 58 cm
Seat depth	43 cm, 46 cm, 49 cm
Back height	40 cm - 47.5 cm
Seat height	Front 40 cm - 51 cm Rear 38 cm - 50 cm
Lower leg length	27 cm - 47 cm