




AERO T

PURE
GEOMETRY.



Take in the curves and lines of TiLite's most advanced aluminum chair. The Aero T's dual-tube frame minimizes weight while maximizing performance. Add to that full adjustability and TiFit, and you've got something special. It's all in the ride.



AERO T

TECHNICAL SPECIFICATIONS.

Frame Features	Aluminum Dual-Tube Rigid TiFit Full Adjustability
Frame Options	Titanium Frame Upgrade SuperLite Edition Armed Forces Edition Power Adaptable Frame Reinforced Frame
Transport Weight*	12 lbs. (5.4 kg)
Weight Limit	265 lbs. (120.2 kg)
Total Chair Width*	23" (58.4 cm)
Seat Width	12" - 20" (30.5 - 51 cm)
Front Seat Width (Seat Taper)	3 Options
Seat Depth	12" - 20" (30.5 - 51 cm)
Front Seat Height	15" - 21" (38 - 53 cm)
Rear Seat Height	13" - 21" (33 - 53 cm)
Camber	0°, 2°, 4°, 6°, 8°, 12°, 15° or 18°
Front Frame Angle	70°, 80°, 85° or 90°
Footrest Width Options	Standard Taper or V-Front End
Center of Gravity Adjustment	5" (12.5 cm)
Front Wheels	23 Options
Rear Wheels	67 Options
Handrims	15 Options
Convenience	Folding Angle Adjustable Seat Back
Aesthetics	5 Aluminum Frame Finishes 19 Frame Colors 6 Color Anodize Packages 7 Ultrasuede Colors 8 Spinergy Wheel Colors
HCPCS Code	K0005
Warranty	Lifetime on Frame One Year on TiLite Components

*Weights / widths based on 16" x 16" (40.5 x 40.5 cm) chair, 0° camber, minimum rear wheel spacing and standard features. Rear wheels not included in transport weight. For reference, two standard 24" (540) rear wheels weigh 7.2 lbs (3.3 kg).



TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. The optimized fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!

To learn more about TiLite options that can take your chair's performance to the next level, visit tilite.com/optimize.



TiLite | The Ultimate Ride™ | © 2014 TiLite, LLC.
800.545.2266 | tilite.com | MS0021_RevA_AeroT_Int
TiLite wheelchairs are shown with available options.