



TRA

ABOVE AND
BEYOND.



Imagine the best, then go beyond, and you will arrive at the TRA. With the power and style of its titanium dual-tube frame, TiFit precision, and full adjustability, the TRA has it all. The result is performance and comfort built around you.



TRA

TECHNICAL SPECIFICATIONS.

Frame Features	Titanium Dual-Tube Rigid TiFit Full Adjustability
Frame Options	TRA Black Edition SuperLite Edition Armed Forces Edition Power Adaptable Frame Reinforced Frame
Transport Weight*	11.1 lbs. (5 kg)
Weight Limit	265 lbs. (120.2 kg)
Total Chair Width*	23" (58.4 cm)
Seat Width	12" - 20" (30.5 - 51 cm)
Front Seat Width (Seat Taper)	3 Options
Seat Depth	12" - 20" (30.5 - 51 cm)
Front Seat Height	15" - 21" (38 - 53 cm)
Rear Seat Height	13" - 21" (33 - 53 cm)
Camber	0°, 2°, 4°, 6°, 8°, 12°, 15° or 18°
Front Frame Angle	70°, 80°, 85° or 90°
Footrest Width Options	Standard Taper or V-Front End
Center of Gravity Adjustment	5" (12.5 cm)
Front Wheels	23 Options
Rear Wheels	67 Options
Handrims	15 Options
Convenience	Folding Angle Adjustable Seat Back
Aesthetics	6 Titanium Frame Finishes 19 Frame Colors 6 Color Anodize Packages 7 Ultrasuede Colors 8 Spinergy Wheel Colors
Warranty	Lifetime on Frame One Year on TiLite Components

*Weights / widths based on 16" x 16" (40.5 x 40.5 cm) chair, 0° camber, minimum rear wheel spacing and standard features. Rear wheels not included in transport weight. For reference, two standard 24" (540) rear wheels weigh 7.2 lbs (3.3 kg).



TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. The optimized fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!

To learn more about TiLite options that can take your chair's performance to the next level, visit tilite.com/optimize.



TiLite | The Ultimate Ride™ | © 2014 TiLite, LLC.
800.545.2266 | tilite.com | MS0019_RevA_TRA_Int
TiLite wheelchairs are shown with available options.