



ZRA
NO
BOUNDARIES.



With full adjustability and a minimalist mono-tube frame, the ZRA dishes out style without boundaries. Add the ZRA's TiFit titanium frame, and you have the complete package. Strength, style, and the ability to set-up and configure your chair just for you.



TECHNICAL SPECIFICATIONS.

Frame Features	Titanium Mono-Tube Rigid TiFit Full Adjustability
Frame Options	ZRA Black Edition SuperLite Edition Armed Forces Edition Heavy Duty Edition Power Adaptable Frame Reinforced Frame
Transport Weight*	10.7 lbs. (4.9 kg)
Weight Limit	265 lbs. (120.2 kg) Standard 300 lbs. (136 kg) Heavy Duty Edition
Total Chair Width*	23" (58.4 cm)
Seat Width	10" - 20" (25 - 51 cm)
Front Seat Width (Seat Taper)	3 Options
Seat Depth	10" - 20" (25 - 51 cm)
Ergonomic Seat	7 Options
Front Seat Height	16" - 21" (40.5 - 53 cm)
Rear Seat Height	12.5" - 20" (32 - 51 cm)
Camber	0°, 2°, 4°, 6°, 8°, 12°, 15° or 18°
Front Frame Angle	70°, 75°, 80°, 85° or 90°
Footrest Width Options	Standard Taper or V-Front End
Center of Gravity Adjustment	6" (15.2 cm)
Front Wheels	23 Options
Rear Wheels	67 Options
Handrims	15 Options
Convenience	Transit Tie-Down Option Folding Angle Adjustable Seat Back Depth Adjustable Back/Frame
Aesthetics	6 Titanium Frame Finishes 19 Frame Colors 6 Color Anodize Packages 7 Ultrasuede Colors 8 Spinergy Wheel Colors
Warranty	Lifetime on Frame One Year on TiLite Components

*Weights / widths based on 16" x 16" (40.5 x 40.5 cm) chair, 0° camber, minimum rear wheel spacing and standard features. Rear wheels not included in transport weight. For reference, two standard 24" (540) rear wheels weigh 7.2 lbs (3.3 kg).



TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. The optimized fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!



TiLite | The Ultimate Ride™ | © 2014 TiLite, LLC.
800.545.2266 | tilite.com | MS0025_RevA_ZRA_Int
TiLite wheelchairs are shown with available options.